

Anti-Ifngr2 Polyclonal Antibody For Wb And Elisa - P38484

Item Number: CM0027036



Introduction

This rabbit anti-IFNGR2 polyclonal antibody is designed for Western blotting (WB) and ELISA applications, with cross-reactivity to human, mouse, and rat, facilitating interferon signaling studies.

[Learn More](#)

Parameter	Description
Product Name	IFNGR2 Rabbit pAb
Remarks/Alias	AF-1; IFGR2; IMD28; IFNGT1; IFNGR2
Species	Human
Gene ID (human)	3460
Gene ID	3460
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 238-337 of human IFNGR2 (NP_005525.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38484
Protein Weight	38kDa
Shipping	Ice bag