

Anti-Ifngr2 Polyclonal Antibody For Wb, Elisa - P38484

Item Number: CM0016716



Introduction

This rabbit polyclonal antibody detects human IFNGR2 (interferon gamma receptor 2), validated for WB and ELISA. Cross-reacts with rat. Ideal for studying interferon gamma signaling and receptor function.

[Learn More](#)

Parameter	Value
Product Name	IFNGR2 Rabbit pAb
Remarks/Alias	AF-1, IFGR2, IMD28, IFNGT1, IFNGR2
Species	Human
Gene ID (Human)	3460
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 221-322 of human IFNGR2 (NP_005525.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P38484
Protein Weight	38kDa
Shipping	Ice bag