

Anti-Ifnl1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8lu54

Item Number: CM0025975



Introduction

Rabbit polyclonal antibody against human IFNL1 (IL-29), validated for WB, IHC-P, ELISA with cross-reactivity to mouse and rat. Ideal for antiviral and immunomodulatory research.

[Learn More](#)

Parameter	Value
Product Name	IFNL1 Rabbit pAb
Remarks/Alias	IL29; IL-29; IFNL1
Species	Human
Gene ID (Human)	282618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human IFNL1 (NP_742152.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IU54
Protein Weight	22kDa
Shipping	Ice bag