

Anti-Ifnl3 Polyclonal Antibody For Ihc-P And Elisa - Q8Izi9

Item Number: CM0015496



Introduction

High-quality rabbit polyclonal antibody targeting human IFNL3 (IL-28B/IL-28C), validated for IHC-P and ELISA, and cross-reactive with mouse and rat. Ideal for investigating type III interferon-mediated antiviral defense in epithelial tissues.

[Learn More](#)

Parameter	Value
Product Name	IFNL3 Rabbit pAb
Remarks/Alias	IL28B; IL28C; IL-28B; IL-28C; IFN-lambda-3; IFN-lambda-4; IFNL3
Species	Human
GeneID (Human)	282617
GeneID	282617
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-196 of human IFNL3 (NP_742151.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZI9
Protein Weight	22kDa
Shipping	Ice bag