

Anti-Ifnl3 Polyclonal Antibody For Wb & Elisa - Q8Izi9

Item Number: CM0027953



Introduction

High-quality rabbit polyclonal antibody against human IFNL3 (IL28B/IL28C). Validated for Western blot and ELISA, with cross-reactivity to human and mouse. For antiviral and immunomodulatory research.

[Learn More](#)

Parameter	Value
Product Name	IFNL3 Rabbit pAb
Remarks/Alias	IL28B; IL28C; IL-28B; IL-28C; IFN-lambda-3; IFN-lambda-4; IFNL3
Species	Human
GeneID (Human)	282617
GeneID	282617
Immunogen	Synthetic peptide within amino acids 1-100 of human IFNL3 (NP_742151.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IZI9
Protein Weight	22kDa
Shipping	Ice bag