

Anti-Ifnlr1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8lu57

Item Number: CM0016960



Introduction

High-quality anti-IFNLR1 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, rat. Suitable for studying antiviral defense and JAK/STAT signaling in epithelial cells.

[Learn More](#)

Parameter	Value
Product Name	IFNLR1 Rabbit pAb
Remarks/Alias	IFNLR; LICR2; IL28RA; CRF2/12; IL-28R1; IFNLR1
Species	Human
GeneID (Human)	163702
GeneID	163702
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 21-228 of human IFNLR1 (NP_734464.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IU57
Protein Weight	58kDa
Shipping	Ice bag