

Anti-Ifrd1 Polyclonal Antibody For Wb, Elisa - 000458

Item Number: CM0018482

Introduction

Rabbit polyclonal antibody targeting human IFRD1 (PC4/TIS7). Validated for Western blot (WB) and ELISA, with cross-reactivity in human, mouse, and rat. Suitable for nerve growth factor (NGF) signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	IFRD1 Rabbit pAb
Remarks/Alias	PC4; TIS7; IFRD1
Species	Human
GeneID (Human)	3475
GeneID	3475
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human IFRD1 (NP_001541.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	50kDa
Shipping	Ice bag