

Anti-Ifrd1 Rabbit Polyclonal Antibody For Wb, Elisa - 000458

Item Number: CM0029723

Introduction

Rabbit polyclonal antibody targeting human IFRD1 (PC4/TIS7) for Western blot and ELISA. Shows cross-reactivity with mouse and rat. Ideal for studying NGF-mediated differentiation and developmental regulation.

[Learn More](#)

Parameter	Value
Product Name	IFRD1 Rabbit pAb
Remarks/Alias	PC4; TIS7; IFRD1
Species	Human
Gene ID (Human)	3475
Gene ID	3475
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00458
Protein Weight	50kDa
Shipping	Ice bag