

Anti-Ift57 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nwb7

Item Number: CM0024721

Introduction

Rabbit polyclonal antibody against human IFT57 (HIPPI, OFD18) validated for WB and ELISA, cross-reactive with mouse and rat. Suitable for cilia and hedgehog signaling research.

[Learn More](#)

Parameter	Value
Product Name	IFT57 Rabbit pAb
Remarks/Alias	HIPPI; OFD18; MHS4R2; ESRRBL1; IFT57
Species	Human
Gene ID	55081
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-429 of human IFT57 (NP_060480.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWB7
Protein Weight	49kDa
Shipping	Ice bag