

Anti-Ift74 Polyclonal Antibody For Wb, Elisa - Q96Lb3

Item Number: CM0017492

Introduction

Intraflagellar transport protein 74 homolog (IFT74) rabbit polyclonal antibody validated for WB and ELISA, reactive with human, mouse, and rat. Ideal for ciliogenesis and flagellogenesis research.

[Learn More](#)

Parameter	Value
Product Name	IFT74 Rabbit pAb
Remarks/Alias	CMG1; BBS22; CCDC2; CMG-1; JBTS40; SPGF58; IFT74
Species	Human
Gene ID (Human)	80173
Gene ID	80173
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-372 of human IFT74 (NP_001092694.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96LB3
Protein Weight	69kDa
Shipping	Ice bag