

Anti-Ift81 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wya0

Item Number: CM0028351

Introduction

Rabbit polyclonal antibody targeting human IFT81 (CDV1), validated for WB and ELISA. Cross-reacts with mouse and rat. Designed for cilia/IFT complex B and SHH signaling research. SWISS Q8WYA0.

[Learn More](#)

Parameter	Value
Product Name	IFT81 Rabbit pAb
Remarks/Alias	DV1; CDV1; CDV-1; CDV1R; CDV-1R; SRTD19; IFT81
Species	Human
GeneID	28981
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-676 of human IFT81 (NP_054774.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WYA0
Protein Weight	80kDa
Shipping	Ice bag