

Anti-Ift88 Monoclonal Antibody For Wb, Elisa - Q13099

Item Number: CM0003869

Introduction

Rabbit monoclonal antibody targeting human IFT88, validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying primary cilium biogenesis and autophagy.

[Learn More](#)

Parameter	Value
Product Name	IFT88 Rabbit mAb
Remarks/Alias	DAF19; TG737; TTC10; hTg737; D13S1056E; IFT88
Species	Human
GeneID (Human)	8100
GeneID	8100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 534-833 of human IFT88. (NP_783195.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13099
Protein Weight	93kDa
Shipping	Ice bag