

Anti-Ift88 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13099

Item Number: CM0013206

Introduction

Rabbit polyclonal IFT88 antibody validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for cilia and autophagy research. (SWISS: Q13099)

[Learn More](#)

Parameter	Value
Product Name	IFT88 Rabbit pAb
Remarks/Alias	DAF19; TG737; TTC10; hTg737; D13S1056E; IFT88
Species	Human
GeneID (Human)	8100
GeneID	8100
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 534-833 of human IFT88 (NP_783195.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13099
Protein Weight	93kDa
Shipping	Ice bag