


Anti-Igf1 Polyclonal Antibody For Wb - P05019

Item Number: CM0020987



Introduction

Rabbit polyclonal antibody against human IGF1 (insulin-like growth factor I) for Western blot. Reacts with mouse and rat. Suitable for studying growth factor signaling, apoptosis, and development.

[Learn More](#)

Parameter	Specification
Product Name	IGF1 Rabbit pAb
Remarks/Alias	IGF; MGF; IGF-I; IGF1
Species	Human
GeneID (Human)	3479
GeneID	3479
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 65-165 of human IGF1 (NP_000609.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Mouse, Rat
SWISS	P05019
Protein Weight	22kDa
Shipping	Ice bag