


Anti-Igf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05019

Item Number: CM0017281



Introduction

Rabbit polyclonal antibody targeting human IGF1 (insulin-like growth factor 1), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse IGF1. Ideal for growth factor signaling, metabolism, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	IGF1 Rabbit pAb
Remarks/Alias	IGF; MGF; IGF1; IGF-I; IGF1
Species	Human
GeneID	3479
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human IGF1 (NP_000609.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05019
Protein Weight	22kDa
Shipping	Ice bag