

Anti-Igf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05019

Item Number: CM0017250



Introduction

High-quality rabbit polyclonal antibody against human IGF1 for Western blot, immunohistochemistry, and ELISA. Cross-reacts with mouse and rat. Applicable in growth factor signaling, metabolism, and neuroscience research.

[Learn More](#)

Product Specifications	Details
Product Name	IGF1 Rabbit pAb
Remarks/Alias	IGF; MGF; IGFI; IGF-I; IGF1
Species	Human
GeneID (human)	3479
GeneID	3479
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human IGF1 (NP_000609.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05019
Protein Weight	22kDa
Shipping	Ice bag