

Anti-Igf1R Rabbit Polyclonal Antibody For Wb, Elisa - P08069

Item Number: CM0018374



Introduction

Rabbit polyclonal antibody targeting human IGF1R (CD221/IGFIR). Validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for studying IGF signaling, cancer biology, and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] IGF1R Rabbit pAb
Remarks/Alias	IGFR; CD221; IGFIR; JTK13; 1R
Species	Human
GeneID	3480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1201-1367 of human IGF1R (NP_000866.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08069
Protein Weight	155kDa
Shipping	Ice bag