


Anti-Igf2 Rabbit Polyclonal Antibody For If/Icc - P01344

Item Number: CM0023642



Introduction

Rabbit polyclonal antibody targeting human IGF2, suitable for IF/ICC. Cross-reacts with mouse and rat. Ideal for studies on fetal growth, glucose metabolism, and tissue differentiation.

[Learn More](#)

Parameter	Specification
Product Name	IGF2 Rabbit pAb
Remarks/Alias	IGF2; C11orf43; GRDF; IGF-II; PP9974; insulin like growth factor 2
Species	Human
GeneID (Human)	3481
GeneID	3481
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 70-180 of human IGF2 (P01344).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC
Cross Reactivity	Human, Mouse, Rat
SWISS	P01344
Protein Weight	20kDa/26kDa
Shipping	Ice bag