

Anti-Igf2Bp1/Imp1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Nzi8

Item Number: CM0017540



Introduction

IGF2BP1/IMP1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat IGF2BP1; ideal for RNA biology, cancer, and neurodevelopment research.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP1/IMP1 Rabbit pAb
Remarks/Alias	IMP1; ZBP1; CRDBP; IMP-1; CRD-BP; VICKZ1; IGF2BP1/IMP1
Species	Human
Gene ID (Human)	10642
Gene ID	10642
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 478-577 of human IGF2BP1/IMP1 (NP_006537.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Nzi8
Protein Weight	63kDa
Shipping	Ice bag