


Anti-Igf2Bp2/Imp2 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Y6M1

Item Number: CM0004970



Introduction

Anti-IGF2BP2/IMP2 Rabbit Monoclonal Antibody for WB, IHC-P, IP, ELISA. Reacts with human IGF2BP2/IMP2 (66/62 kDa), an RNA-binding protein regulating mRNA stability and transport.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP2/IMP2 Rabbit mAb
Remarks/Alias	IMP2; IMP-2; VICKZ2; IGF2BP2/IMP2
Species	Human
GeneID (Human)	10644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-595 of human IGF2BP2/IMP2 (NP_001007226.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human
SWISS	Q9Y6M1
Protein Weight	66 kDa/62 kDa
Shipping	Ice bag