


Anti-Igf2Bp2/Imp2 Polyclonal Antibody For If/Icc, Elisa - Q9Y6M1

Item Number: CM0013165



Introduction

This rabbit polyclonal antibody targets human IGF2BP2/IMP2 (also known as IMP-2, VICKZ2), validated for IF/ICC and ELISA. Cross-reacts with mouse and rat. Study its function in mRNA transport and translational control.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP2/IMP2 Rabbit pAb
Remarks/Alias	IMP2; IMP-2; VICKZ2; IGF2BP2/IMP2
Species	Human
GeneID	10644 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 500-599 of human IGF2BP2/IMP2 (NP_006539.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6M1
Protein Weight	66kDa
Shipping	Ice bag