


Anti-Igf2Bp2/Imp2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y6M1

Item Number: CM0026239



Introduction

High-quality rabbit polyclonal antibody against IGF2BP2/IMP2, a key RNA-binding protein involved in mRNA transport and stability. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human and mouse. Ideal for cancer and developmental research.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP2/IMP2 Rabbit pAb
Remarks/Alias	IMP2; IMP-2; VICKZ2; IGF2BP2/IMP2
Species	Human
GeneID (human)	10644
GeneID	10644
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 500-599 of human IGF2BP2/IMP2 (NP_006539.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6M1
Protein Weight	66kDa
Shipping	Ice bag