


Anti-Igf2Bp2/Imp2 Rabbit Polyclonal Antibody For Wb, If, Elisa - Q9Y6M1

Item Number: CM0022859



Introduction

Rabbit polyclonal antibody against human IGF2BP2/IMP2. Validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Targets a ~66kDa protein involved in mRNA stability and transport.

[Learn More](#)

Parameter	Description
Product Name	IGF2BP2/IMP2 Rabbit pAb
Remarks/Alias	IMP2; IMP-2; VICKZ2; IGF2BP2/IMP2
Species	Human
GeneID	10644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-599 of human IGF2BP2/IMP2 (NP_006539.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6M1
Protein Weight	66kDa
Shipping	Ice bag