


Anti-Igf2Bp3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O00425

Item Number: CM0029972



Introduction

Rabbit polyclonal antibody against human IGF2BP3, validated for WB, IHC-P, IF/ICC, IP, ELISA with cross-reactivity to mouse and rat. Ideal for RNA-binding protein research in cancer and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP3 Rabbit pAb
Remarks/Alias	KOC, CT98, IMP3, KOC1, IMP-3, VICKZ3, IGF2BP3
Species	Human
GeneID	10643 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-579 of human IGF2BP3 (NP_006538.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00425
Protein Weight	64kDa
Shipping	Ice bag