


Anti-Igf2Bp3/Imp3 Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - O00425

Item Number: CM0006563



Introduction

Recombinant rabbit monoclonal antibody against human IGF2BP3/IMP3, validated for WB, IP, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for cancer and RNA-binding protein research.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP3/IMP3 Rabbit mAb
Remarks/Alias	KOC; CT98; IMP3; KOC1; IMP-3; VICKZ3; IGF2BP3/IMP3
Species	Human
GeneID (Human)	10643
GeneID	10643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-579 of human IGF2BP3/IMP3 (NP_006538.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00425
Protein Weight	64 kDa
Shipping	Ice bag