


Anti-Igfals Polyclonal Antibody For Wb, Elisa - P35858

Item Number: CM0013889



Introduction

Rabbit anti-human IGFALS polyclonal antibody validated for WB and ELISA, with cross-reactivity to mouse and rat. IGFALS, a secreted plasma protein, regulates IGF bioavailability via ternary complex formation. Ideal for studying protein interactions and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	IGFALS Rabbit pAb
Remarks/Alias	ALS; ACLSD; IGFALS
Species	Human
GeneID (Human)	3483
GeneID	3483
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 396-605 of human IGFALS (NP_004961.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35858
Protein Weight	66kDa
Shipping	Ice bag