

Anti-Igfbp1 Monoclonal Antibody For Wb, Elisa - P08833

Item Number: CM0002865



Introduction

Rabbit monoclonal antibody against human IGFBP1, suitable for Western blot and ELISA. IGFBP1 regulates IGF bioavailability and cell migration; useful for growth factor research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	IGFBP1 Rabbit mAb
Remarks/Alias	AFBP, IBP1, PP12, IGF-BP25, hIGFBP-1, IGFBP1
Species	Human
GeneID (Human)	3484
GeneID	3484
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 160-259 of human IGFBP1 (P08833).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08833
Protein Weight	28kDa
Shipping	Ice bag