

Anti-Igfbp1 Rabbit Polyclonal Antibody For Wb, Elisa - P08833

Item Number: CM0017045



Introduction

IGFBP1 Rabbit pAb is a polyclonal antibody targeting human insulin-like growth factor-binding protein 1, validated for Western blot and ELISA. Suitable for research on IGF signaling, cell migration, and apoptosis regulation.

[Learn More](#)

Parameter	Value
Product Name	IGFBP1 Rabbit pAb
Remarks/Alias	AFBP; IBP1; PP12; IGF-BP25; hIGFBP-1; IGFBP1
Species	Human
GeneID (Human)	3484
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 150-250 of human IGFBP1 (NP_000587.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P08833
Protein Weight	28kDa
Shipping	Ice bag