


Anti-Igfbp2 Monoclonal Antibody For Wb, Ihc-P, Elisa - P18065

Item Number: CM0009546



Introduction

Rabbit monoclonal antibody targeting human IGFBP2. Validated for WB, IHC-P, and ELISA. Useful for research on insulin-like growth factor binding, cell migration, and mitochondrial apoptosis.

[Learn More](#)

Parameter	Value
Product Name	IGFBP2 Rabbit mAb
Remarks/Alias	IBP2; IGF-BP53; IGFBP2
Species	Human
Gene ID	3485
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 226-325 of human IGFBP2 (P18065)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P18065
Protein Weight	35kDa
Shipping	Ice bag