


## Anti-Igfbp2 Rabbit Polyclonal Antibody For Wb And Elisa - P18065

Item Number: CM0018020



### Introduction

Rabbit polyclonal antibody targeting human IGFBP2, for WB and ELISA. Cross-reacts with mouse and rat. Recognizes a ~35 kDa protein. Ideal for investigating IGF-mediated proliferation, differentiation, and migration.

[Learn More](#)

| Parameter        | Value                                                                                                  |
|------------------|--------------------------------------------------------------------------------------------------------|
| Product Name     | IGFBP2 Rabbit pAb                                                                                      |
| Remarks/Alias    | IBP2; IGF-BP53; IGFBP2                                                                                 |
| Species          | Human                                                                                                  |
| Gene ID (Human)  | 3485                                                                                                   |
| Gene ID          | 3485                                                                                                   |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 100-200 of human IGFBP2 (NP_000588.2) |
| Source           | Rabbit                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                  |
| Application      | WB, ELISA                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                      |
| SWISS            | P18065                                                                                                 |
| Protein Weight   | 35kDa                                                                                                  |
| Shipping         | Ice bag                                                                                                |