

Anti-Igfbp4 Rabbit Polyclonal Antibody For Wb, Elisa - P22692

Item Number: CM0015630



Introduction

IGFBP4 Rabbit Polyclonal Antibody validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying IGF bioavailability and growth regulation.

[Learn More](#)

Parameter	Value
Product Name	IGFBP4 Rabbit pAb
Remarks/Alias	BP-4; IBP4; IGFBP-4; HT29-IGFBP; IGFBP4
Species	Human
GeneID (Human)	3487
GeneID	3487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 178-258 of human IGFBP4 (NP_001543.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P22692
Protein Weight	28kDa
Shipping	Ice bag