

Anti-Igfbp4 Rabbit Polyclonal Antibody - P22692

Item Number: CM0020386



Introduction

Rabbit polyclonal antibody against human IGFBP4, validated for WB, IHC-P, IF/ICC, and ELISA. Detects IGFBP4 in human, mouse, and rat samples. Suitable for studying IGF signaling and growth regulation.

[Learn More](#)

Parameter	Value
Product Name	IGFBP4 Rabbit pAb
Remarks/Alias	BP-4; IBP4; IGFBP-4; HT29-IGFBP; IGFBP4
Species	Human
Gene ID	3487
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-258 of human IGFBP4 (NP_001543.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22692
Protein Weight	28kDa
Shipping	Ice bag