

Anti-Igfbp7 Rabbit Polyclonal Antibody For Wb, Elisa - Q16270

Item Number: CM0016916



Introduction

Rabbit polyclonal antibody against human IGFBP7, validated for WB and ELISA. Recognizes mouse and rat IGFBP7. Ideal for studying insulin-like growth factor binding, cell adhesion, and angiogenesis.

[Learn More](#)

Parameter	Specification
Product Name	IGFBP7 Rabbit pAb
Remarks/Alias	AGM; PSF; TAF; FSTL2; IBP-7; MAC25; IGFBP-7; RAMSVPS; IGFBP-7v; IGFBPRP1; IGFBP7
Species	Human
GeneID (human)	3490
GeneID	3490
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-282 of human IGFBP7 (NP_001544.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16270
Protein Weight	29kDa
Shipping	Ice bag