

Anti-Igfbpl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wx77

Item Number: CM0023200



Introduction

Rabbit polyclonal antibody targeting human IGFBPL1 (IGFBP-rP4). Validated in WB and ELISA; cross-reacts with mouse and rat. Suitable for studying IGF signaling and tumor suppression. 29kDa.

[Learn More](#)

Parameter	Value
Product Name	IGFBPL1 Rabbit pAb
Remarks/Alias	IGFBPRP4; IGFBP-RP4; bA113O24.1; IGFBPL1
Species	Human
Gene ID	347252
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 150-250 of human IGFBPL1 (NP_001007564.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WX77
Protein Weight	29kDa
Shipping	Ice bag