


Anti-Ighd Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P01880

Item Number: CM0002304



Introduction

Anti-IGHD Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA. Targets human immunoglobulin heavy constant delta chain (P01880) for B-cell and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	IGHD Rabbit mAb
Species	Human
GeneID (human)	3495
GeneID	3495
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human IGHG (P01880).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P01880
Protein Weight	42kDa
Shipping	Ice bag