


Anti-Ighmbp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P38935

Item Number: CM0020902



Introduction

Rabbit polyclonal antibody against human IGHMBP2, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for research on RNA helicase, ribosome biogenesis, and neurodegenerative diseases.

[Learn More](#)

Parameter	Value
Product Name	IGHMBP2 Rabbit pAb
Remarks/Alias	HCSA; HMN6; CATF1; CMT2S; SMARD1; SMUBP2; ZFAND7; IGHMBP2
Species	Human
Gene ID (Human)	3508
Gene ID	3508
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-720 of human IGHMBP2 (NP_002171.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38935
Protein Weight	109kDa
Shipping	Ice bag