

Anti-Igsf3 Rabbit Polyclonal Antibody For Wb, Elisa - 075054

Item Number: CM0016289

Introduction

This rabbit polyclonal antibody against human IGSF3 (EWI-3) is validated for Western blotting and ELISA. Suitable for studying immunoglobulin superfamily member 3 expression in placenta, kidney, and lung tissues.

[Learn More](#)

Parameter	Value
Product Name	IGSF3 Rabbit pAb
Remarks/Alias	V8; LCDD; EWI-3; IGSF3
Species	Human
Gene ID	3321
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-562 of human IGSF3 (NP_001007238.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75054
Protein Weight	135kDa
Shipping	Ice bag