

Anti-Ihh Rabbit Polyclonal Antibody For Wb And Elisa - Q14623

Item Number: CM0024884



Introduction

Rabbit polyclonal antibody targeting IHH (Indian hedgehog protein), suitable for WB and ELISA. Reacts with human, mouse, and rat samples. Useful for research on bone and cartilage development, cell signaling, and related fields.

[Learn More](#)

Parameter	Value
Product Name	IHH Rabbit pAb
Remarks/Alias	BDA1; HHG2; IHH
Species	Human
Gene ID	3549
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 202-411 of human IHH (NP_002172.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14623
Protein Weight	45kDa
Shipping	Ice bag