


Anti-Ik Polyclonal Antibody For Wb, Elisa - Q13123

Item Number: CM0014967



Introduction

Rabbit polyclonal antibody targeting human IK (Red/RER). Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for spliceosome and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	IK Rabbit pAb
Remarks/Alias	RED; RER; CSA2; IK
Species	Human
Gene ID (Human)	3550
Gene ID	3550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human IK (NP_006074.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13123
Protein Weight	66kDa
Shipping	Ice bag