


Anti-Ikaros Rabbit Polyclonal Antibody For Wb, Elisa - Q13422

Item Number: CM0024698



Introduction

Anti-Ikaros Rabbit Polyclonal Antibody (Q13422) validated for WB and ELISA in human. Targets the Ikaros transcription factor involved in lymphocyte and erythroid differentiation. Ideal for immunodetection studies.

[Learn More](#)

Parameter	Value
Product Name	Ikaros Rabbit pAb
Remarks/Alias	IK1; LYF1; LyF-1; CVID13; IKAROS; PPP1R92; PRO0758; ZNFN1A1; Hs.54452; Ikaros
Species	Human
GeneID (Human)	10320
GeneID	10320
Immunogen	Recombinant protein A synthetic peptide corresponding to a sequence within amino acids 296-457 of human Ikaros(NP_006051.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13422
Protein Weight	58kDa
Shipping	Ice bag