

Anti-Ikaros Rabbit Polyclonal Antibody For Wb - Q13422

Item Number: CM0022475



Introduction

Rabbit polyclonal antibody against human Ikaros (IKZF1), validated for WB, IF/ICC, ELISA, ChIP. Cross-reacts with mouse and rat. Ideal for hematopoietic differentiation and lymphocyte development research.

[Learn More](#)

Parameter	Value
Product Name	Ikaros Rabbit pAb
Remarks/Alias	IK1; LYF1; LyF-1; CVID13; IKAROS; PPP1R92; PRO0758; ZNFN1A1; Hs.54452; Ikaros
Species	Human
GeneID (human)	10320
GeneID	10320
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-270 of human Ikaros (NP_001207694.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13422
Protein Weight	58kDa
Shipping	Ice bag