

Anti-Ikba Polyclonal Antibody For Wb, Elisa - P25963

Item Number: CM0013236



Introduction

High-quality IkBα Rabbit pAb validated for Western blot and ELISA. Detects endogenous NF-kappa-B inhibitor alpha in human, mouse, and rat samples. Ideal for NF-κB pathway and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IkBα Rabbit pAb
Remarks/Alias	IKBA; MAD-3; NFKBI; EDAID2; IkBα
Species	Human
Gene ID (Human)	4792
Gene ID	4792
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 166-265 of human IkBα (NP_065390.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25963
Protein Weight	36kDa
Shipping	Ice bag