

Anti-Ikba Polyclonal Antibody For Wb, Ip, Elisa - P25963

Item Number: CM0016911



Introduction

Rabbit polyclonal antibody against human IκBα (NF-kappa-B inhibitor alpha), suitable for WB, IP, ELISA. Cross-reacts with mouse and rat. Ideal for studying NF-κB signaling, inflammation, and immune responses.

[Learn More](#)

Parameter	Value
Product Name	IκBα Rabbit pAb
Remarks/Alias	IKBA; MAD-3; NFKBI; EDAID2; IκBα
Species	Human
GeneID (Human)	4792
GeneID	4792
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 201-300 of human IκBα (NP_065390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25963
Protein Weight	36kDa
Shipping	Ice bag