

Anti-Ikba Rabbit Monoclonal Antibody For Wb, Elisa - Q63746

Item Number: CM0011244

Introduction

Rabbit monoclonal antibody targeting IkBα (NFKBIA). Validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Suitable for NF-kappa-B signaling studies. Predicted molecular weight: 35 kDa.

[Learn More](#)

Parameter	Value
Product Name	IkBα Rabbit mAb
Remarks/Alias	RL/IF-1; IkBα
Species	Rat
GeneID (Human)	4792
GeneID	25493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-314 of Rat IkBα (NP_001099190.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q63746
Protein Weight	35 kDa
Shipping	Ice bag