

Anti-Ikbip Polyclonal Antibody For Wb, If/Icc, Elisa - Q70Uq0

Item Number: CM0028787

Introduction

Rabbit polyclonal antibody targeting IKBIP (inhibitor of nuclear factor kappa-B kinase-interacting protein), a p53-regulated pro-apoptotic factor. Validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	IKBIP Rabbit pAb
Remarks/Alias	IKIP; IKBIP
Species	Human
GeneID (human)	121457
GeneID	121457
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 115-350 of human IKBIP (NP_963906.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q70UQ0
Protein Weight	39kDa
Shipping	Ice bag