


Anti-Ikbkap Polyclonal Antibody For Wb, Ihc-P, Elisa - 095163

Item Number: CM0018009



Introduction

High-quality anti-IKBKAP rabbit polyclonal antibody validated for WB, IHC-P, and ELISA with human, mouse, and rat reactivity. Ideal for studying elongator complex, tRNA modification, and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	IKBKAP Rabbit pAb
Remarks/Alias	FD; DYS; IKAP; IKI3; TOT1; IKBKAP
Species	Human
GeneID (Human)	8518
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1093-1332 of human IKBKAP (NP_003631.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95163
Protein Weight	150kDa
Shipping	Ice bag