


Anti-Ikk Beta Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Elisa - O14920

Item Number: CM0008520



Introduction

High-specificity recombinant rabbit monoclonal antibody against human IKK beta. Validated for Western blot, immunoprecipitation, IF/ICC, and ELISA. Ideal for NF- κ B signaling and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IKK beta Rabbit mAb
Remarks/Alias	IKK2, IKKB, IKK-2, IMD15, IMD15A, IMD15B, NFKBIKB, IKK-beta
Species	Human
GeneID	3551
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 261-411 of human IKK beta (NP_001547.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O14920
Protein Weight	29-87 kDa
Shipping	Ice bag