


Anti-Ikk α Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O15111

Item Number: CM0022893



Introduction

Rabbit polyclonal antibody against human IKK α , detecting an ~85kDa band. Validated for WB, IHC-P, IF/ICC, IP, ELISA with human, mouse, rat samples. Ideal for studying canonical and non-canonical NF- κ B pathways.

[Learn More](#)

Parameter	Value
Product Name	IKK α Rabbit pAb
Remarks/Alias	BPS2; IKK1; IKKA; IKBKA; TCF16; NFKBKA; IKK-alpha; IKK α
Species	Human
GeneID (Human)	1147
GeneID	1147
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human IKK α (NP_001269.3).
Source	Rabbit
Category	Polyclonal
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15111
Protein Weight	85kDa
Shipping	Ice bag