

Anti-Ikkε (Ikbke) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14164

Item Number: CM0013026

Introduction

Rabbit polyclonal antibody targeting human IKKε (IKBKE), suitable for Western blot, immunohistochemistry, and ELISA. Detects human, mouse, and rat IKKε. An essential tool for studying innate immunity and inflammatory signaling.

[Learn More](#)

Parameter	Specification
Product Name	IKKε Rabbit pAb
Remarks/Alias	IKKE; IKKI; IKK-E; IKK-i; IKKε
Species	Human
GeneID (human)	9641
GeneID	9641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 495-716 of human IKKε (NP_054721.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14164
Protein Weight	80kDa
Shipping	Ice bag