

Anti-Ikkε (Ikbke) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q14164

Item Number: CM0018577

Introduction

IKKε (IKBKE) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects endogenous IKKε in human, mouse, and rat samples. Suitable for studying antiviral immunity and inflammatory signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	IKKε Rabbit pAb
Remarks/Alias	IKKE; IKKI; IKK-E; IKK-i; IKKε
Species	Human
GeneID (Human)	9641
GeneID	9641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 497-716 of human IKKε (NP_054721.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14164
Protein Weight	80kDa
Shipping	Ice bag